

杭州桐庐康尔医疗器械有限公司

HANGZHOU TONGLU KANGER MEDICAL INSTRUMENT CO.,LTD

杭州桐庐富尔克医疗器械有限公司

HANGZHOU TONGLU FUERKE MEDICAL INSTRUMENT CO.,LTD

通过ISO9001、ISO13485国际质量管理体系认证
通过欧盟CE认证



杭州桐庐康尔医疗器械有限公司

HANGZHOU TONGLU KANGER MEDICAL INSTRUMENT CO.,LTD

杭州桐庐富尔克医疗器械有限公司

HANGZHOU TONGLU FUERKE MEDICAL INSTRUMENT CO.,LTD

营销部电话：0571-64250383

地址：浙江省杭州市桐庐县县城阆苑路108号

电话：0571-64266318 64266317 64266316

传真：0571-64266318

邮编：311500

<http://www.hz-keqx.com>

<http://hz-keqx.en.alibaba.com>

E-mail: web@hz-keqx.com

Address: 108 Langyuan road, Tonglu Economic Development Zone, Zhejiang province.

Tel: +86571-64266317

Tel of sales Department: +86 571-64250383

Fax: +86 571-64266318

Zip code: 311500



公司网站



Website in English version

19-20中英文版

科学 诚信 创新 分享

Science Integrity Innovation Share

杭州桐庐康尔医疗器械有限公司
HANGZHOU TONGLU KANGER MEDICAL INSTRUMENT CO.,LTD
杭州桐庐富尔克医疗器械有限公司
HANGZHOU TONGLU FUERKE MEDICAL INSTRUMENT CO.,LTD

企业简介

杭州桐庐康尔医疗器械有限公司是一家浙江省高新技术企业,创建于1999年4月,位于美丽的富春江畔省级经济开发区—桐庐经济开发区内,厂区占地面积6500平方米,是一家集开发、生产、销售于一体,专业生产腹腔镜、胸腔镜、妇科腹腔镜等手术器械的厂家。

公司产品品种多样、规格齐全、质量可靠、款式新颖,售后服务优良,产品畅销全国各大医院。经临床使用,得到专家、教授的一致好评。本公司本着“科学、诚信、创新、分享”的经营理念,以市场需求为导向,以优质的产品为目标,优秀的人才为基础,优良的服务为后盾,开拓创新,不断进取。


桐庐富尔克医疗器械有限公司在康尔公司的基础上成立的,公司本着“质量为本,服务满意”的宗旨,主要销售医用内窥镜和配套用手术器械以及图像显示仪等产品。欢迎广大医疗界人士来函交流指导与洽谈业务。



Company Profile

Hangzhou Tonglu Kanger medical instrument Co.,Ltd has been designated as one of the high-tech enterprise in Zhejiang Province. It was founded in April, 1999, located at a Provincial Economic development zone of Tonglu. It's along the Fuchun River, covering 6500 square meters of the factory area. The company specializing in developing, manufacturing and selling surgical instruments of Laparoscope, Thoracoscope, Gynecological Laparoscope and so on .

Our products have the variety specifications, update models, providing superior quality with the high standard services at competitive price. We are in line with "science, good faith, innovation, sharing" as the company's Managing motto. Putting the market demands as the guidance, the high quality product as the goal, the professional elitists as a foundation, the fine service as the backing, to development innovation and unceasingly motivated.



目录 List

一. 腹腔镜手术器械 Laparoscope instruments	
手术器械拆卸图	Instrument to unload 3
腹腔镜内窥镜	Laparoscope 3
松夹器	Loose clip forces 4
钛夹钳、钛夹、施夹钳	Clip applicator, TI clip, Hem-o-lok 5
施夹器、连接器、导光束	Hem-o-lok, Endoscopy Interface, Light transmitting bundle 6
双极电凝钳	Bipolar forceps 7
弹簧抓钳、荷包钳、档板、牵引	Spring grasping forceps, Auto suture purstring, Baffle, Pull 9
金手指拉钩	Goldfinger retractor 10
持针钳	Needle-holding forceps 11
电凝器械、	Electrode instruments 12
冲洗吸引器、胆管刀、穿刺器	Suction & Irrigation tube, Bile duct knife, Bile duct needle 13
其它手术器械	Other instruments 14
撑开器	Spreader 14
KE-104胸腔手术器械	Thoracic operation equipment 15
剪刀、分离钳	Scissors, Dissecting forceps 19
其它抓钳	Other forceps 20
其它抓钳	Other forceps 21
其它抓钳	Other forceps 23
一次性手术器械	Disposable surgical instrument 24
穿刺器、转换器	Trocar, Reducer 25
KE-102微型腹腔镜手术器械	Miniature instruments for Laparoscopic surgery 26
KE-103小切口手术器械	Instruments for small incision operation 27
KE-104胸腔镜手术器械	Instruments for Thoracoscopic surgery 28
KE-106腹腔镜甲状腺手术器械	Thyroid operating instruments 28
KE-107腹腔镜胆道取石手术器械	Instruments of Biliary tract fetch 28
KE-101DK单孔(经脐部)腹腔镜手术器械	Single-hole(After umbilical) 29
二. 妇科手术器械 Gynaecology	
KE-201A电动子宫切除器	Electric uterus cutter 30
KE-201C妇科腹腔镜手术器械(选购件)	Instruments of Gynaecology Laparoscope 31
KE-201C妇科腹腔镜手术器械(选购件)	Instruments of Gynaecology Laparoscope 32
KE-201C妇科腹腔镜手术器械(选购件)	Instruments of Gynaecology Laparoscope 32
KE-201B子宫切除器械(阴式)	Uterus-cutter(process vagina) 33
KE-202宫腔镜	Hysteroscope 33
KE-203宫腔镜(检查镜)	Cavitas inspectoscope 33
KE-204宫腔电切镜	Hystero Resectoscope 34
三. 泌尿科手术器械 Instruments for Urological Department	
KE-301尿道膀胱镜	Cysto-Urethroscope 34
KE-304膀胱碎石镜	Lithotriposcope 35
KE-305尿道切开镜	Urethrotome telescope 35
四. 耳鼻喉科手术器械 Instruments for E.N.T Department	
KE-501鼻窦镜	Sinoscope 36
KE-503鼻窦电动手术刀	Nose electric drill 36
KE-502支撑喉镜	Hold Laryngoscope 37
KE-506食道镜	Esophagoscope instruments 37
KE-505耳科手术器械包	Otoscope instruments 37
KE-504小儿支气管镜	Bronchoscope 38
五. 骨科手术器械 Instruments for Orthopaedics	
KE-401关节镜	Arthroscopy instruments 38
KE-402椎间盘镜手术器械	Endoscopic Discectomy system 38
六. 其它附属设备 Other equipments	
KE-601数字内窥镜图像显示系统	Endoscope display systems 39
KE-610耳鼻喉科诊疗台系列	E.N.T Therapeutic equipments 39
KE-611模拟操作箱	Simulated training box 39
KE-603医用加压器(冲吸泵)	Suction & Irrigation apparatus 40
KE-604医用加压器(灌注泵)	Pump 40
KE-605冷光源	Cold light source 41
KE-606器械清洗成套设备	Cleaning device for instruments 41
KE-607气腹机	Insufflator 41
KE-608高频电刀	High frequency electric-cutter 42
内窥镜消毒盒	Sterilising trays 42


手术钳拆卸图 Instrument to unload




腹腔镜内窥镜 Laparoscopes




施夹器 Hem-o-lok




松夹器 Loose clip forceps



钛夹钳. 钛夹 Clip applicator. TI Clip




施夹器 (开放式) Hem-o-lok




施夹器 Hem-o-lok




连接器 Endoscopy Interface 导光束 Light transmitting bundle



双极电凝钳 Bipolar forceps



双极电凝钳 Bipolar forceps



弹簧抓钳 Spring grasping forceps



荷包钳 Auto suture purstring



挡板、牵引 Baffle、Pull





金手指拉钩 Goldfinger retractor




持针钳 Needle-holding Forceps



持针钳 Needle-holding Forceps




电凝器械 Electrode instruments




冲洗吸引器 Section & Irrigation tube



胆管刀 Bile duct knife 穿刺针 Bile Duct Needle



其它手术器械 Other instruments



撑开器 Spreader



KE-104胸腔手术器械 Thoracic operation equipment




104.101 双关节分离钳



104.101 A 双关节分离钳 (弧高15)




104.101 B 双关节分离钳 (弧高20)



104.101 C 双关节分离钳 (弧高25)



104.101 D 双关节分离钳 (弧高25,直角)



104.101 E 双关节分离钳

注：双关节滑杆式分离钳（又名止血钳），可以柔和、安全夹住血管，头部设计有利于安全有效地分离血管或支气管，手柄分带锁扣和不带锁扣两种



104.102 双关节持针钳



104.102 A 双关节持针钳 (直头)



104.102 B 双关节持针钳 (弯头)

注：碳化钨合金镶嵌钳口非常安全有效牢固地加持缝针，直行和弯行钳口能根据医生的使用习惯选择



104.103 双关节剪刀



104.103 A 双关节剪刀 (直头)



104.103 B 双关节剪刀 (弯头)

注：双关节滑杆式剪刀制造一个干净的切口，圆钝的头部段防损伤到组织。这些剪刀是由经验丰富的技工手工制作以确保卓越的平衡感及触感，可保证与控制精准的剪切精度。微弯行的剪刀分弯头和直头二种能根据医生需要选择

KE-104胸腔手术器械 Thoracic operation equipment



104.104 双关节卵圆钳



104.104 A 双关节卵圆钳 (8mm有齿)



104.104 B 双关节卵圆钳 (10mm有齿)



104.104 C 双关节卵圆钳 (12mm有齿)



104.104 D 双关节卵圆钳 (8mm无齿)



104.104 E 双关节卵圆钳 (10mm无齿)



104.104 F 双关节卵圆钳 (12mm无齿)

注：双关节滑杆式卵圆钳安全无损地抓起肺部组织，不会造成撕裂或损伤。椭圆形的钳口能够在分离、解剖淋巴结时安全地夹持淋巴结，卵圆形的椭圆形环状钳口能够将折叠的海绵作为海绵垫抓起，用于钝性分离也可根据不同需要、不同病人、不同尺寸与不同组织结构可进行选择



104.105 双关节肺叶钳


注：双关节肺叶钳使用功能类同与双关节滑杆式卵圆钳




104.106 双关节夹持钳

注：双关节滑杆式夹持钳（又名镊子）用于缝合组织时传递缝针，抽拉缝线。方便在针穿过组织时将它夹住，并将针递给持针器。

KE-104胸腔手术器械 Thoracic operation equipment



KE-104胸腔手术器械 Thoracic operation equipment



剪刀、分离钳 Scissors Dissecting forceps



101.006C胆囊抓钳
Grasping forceps Φ10X330



101.013B弯分离钳圆头 (17mm)
Dissecting forceps(curved)
Φ5X330



101.013K弯分离钳横齿 (12mm)
Dissecting forceps(curved)
Φ5X330



101.010大型分离钳
Dissecting forceps(large)
Φ10X330



101.013C弯分离钳方头 (17mm)
Dissecting forceps(curved)
Φ5X330



101.013L弯分离钳横齿
Dissecting forceps(curved)
Φ5X330



101.011A直角分离钳
Dissecting forceps(90°)
Φ10X330



101.013D弯分离钳 (14mm)
Dissecting forceps(curved)
Φ5X330



101.014大型弯剪刀
Scissors(large curved)
Φ10X330



101.011B直角分离钳
Dissecting forceps(90°)
Φ5X330



101.013E弯分离钳钝头形 (16mm)
Dissecting forceps(curved)
Φ5X330



101.015大型弯剪刀
Scissors(large curved)
Φ10X330



101.012直分离钳
Dissecting forceps(straight)
Φ5X330



101.013F弯分离钳圆头 (20mm)
Dissecting forceps(curved)
Φ5X330



101.016单开直剪刀
Scissors(single-action straight)
Φ5X330



101.012B尖头钳
Probe forceps
Φ5X330



101.013G弯分离钳花齿 (20mm)
Dissecting forceps(curved)
Φ5X330



101.017双开直剪刀
Scissors(double-action straight)
Φ5X330



101.013A弯分离钳 (12mm)
Dissecting forceps(curved)
Φ5X330



101.013H弯分离钳圆头 (25mm)
Dissecting forceps(curved)
Φ5X330



101.018单开弯剪刀
Scissors(single-action curved)
Φ5X330

其它抓钳 Other forceps



101.019A双开弯剪刀
Scissors(double-action curved)
Φ5X330



101.036B鸭嘴抓钳(粗齿碗形)
Ischaemum grasping forceps
Φ5X330



101.040D肠钳
Bowel forceps
Φ5X330



101.019B双开弯剪刀(长)
Scissors(double-action curved)
Φ5X330



101.037A中空横齿抓钳
Mid-hollow grasping forceps
Φ5X330



101.041A粗齿无损伤抓钳
Tumour grasper with teeth
Φ5X330



101.022钩剪刀
Hook scissors
Φ5X330



101.037B中空竖齿抓钳
Mid-hollow grasping forceps
Φ5X330



101.041B粗齿抓钳(九齿)
Tumour grasper with teeth
Φ5X330



101.035胆囊抓钳
Gallbladder grasping forceps
Φ5X330



101.040A肠钳(粗齿)
Bowel forceps
Φ5X330



101.041C粗齿抓钳(六孔)
Tumour grasper with teeth
Φ5X330



101.035B胆囊抓钳(三角头)
Gallbladder grasping forceps
Φ5X330



101.040B肠钳(细齿)
Bowel forceps
Φ5X330



101.042A输卵管抓钳
Oviduct forceps
Φ5X330



101.036鸭嘴抓钳
Duck-jaw forceps
Φ5X330



101.040B2肠钳(细齿中空)
Bowel forceps
Φ5X330



101.042B输卵管抓钳
Oviduct forceps
Φ5X330



101.036A鸭嘴抓钳(长头)
Ischaemum grasping forceps
Φ5X330



101.040C肠钳
Bowel forceps
Φ5X330



101.042C输卵管抓钳
Oviduct forceps
Φ5X330

其它抓钳 Other forceps



101.042D V输卵管抓钳
Oviduct forceps
Φ5X330



101.042F 输卵管抓钳(O型)
Oviduct forceps
Φ5X330



101.043A 双动鼠齿抓钳(粗齿)
Claw forceps
Φ5X330



101.043B1 双动鼠齿抓钳(4+2细齿)
Claw forceps
Φ5X330



101.043B2 双动鼠齿抓钳(3+2)
Claw forceps
Φ5X330



101.043C 双动细齿鼠牙(2+2)
Claw forceps
Φ5X330



101.043D 单动细齿鼠牙(4+2)
Claw forceps
Φ5X330




101.044 钝头抓钳
Grasping forceps(blunt)
Φ5X330




101.044B 钝头抓钳(粗齿)
Grasping forceps(blunt)
Φ5X330




101.044A 钝头抓钳(有齿)
Appendiceal forceps
Φ5X330



101.044B 钝头抓钳(粗齿)
Grasping forceps(blunt)
Φ5X330



101.045 胆道造影抓钳
Cholangiography forceps
Φ5X330



101.046 重型抓钳
Max atraumatic forceps
Φ5X330



101.047A 倒齿抓钳
Rerersing tooth forceps
Φ5X330



101.047B 倒齿抓钳(短头)
Rerersing tooth forceps
Φ5X330



101.047C 倒齿抓钳
Rerersing tooth forceps
Φ5X330



101.048 阑尾抓钳
Appendiceal forceps
Φ5X330




101.048A 阑尾抓钳(有齿)
Appendiceal forceps
Φ5X330



101.052 S型打结钳
S-knot forceps
Φ5X330



101.053 B型打结钳
B-knot forceps
Φ5X330



101.054 弯分离型打结钳
Dissecting-knot forceps(curved)
Φ5X330



101.064 带针活检钳
Biopsy forceps(with needle)
Φ5X330



101.065 活检钳
Biopsy forceps
Φ5X330



101.066B 三叶扇形钳
Fan retractor forceps
Φ5X330

其它抓钳 Other forceps



101.066B 三叶扇形钳
Fan retractor forceps
Φ5X330



101.067C 扇形钳
Sector clamp
Φ5X330



101.076B 肺叶钳
Lung lobe forceps
Φ5X330



101.077 海绵钳
Absorbent cotton forceps
Φ5X330



101.078 系膜钳
Mesentery forceps
Φ5X330



101.078B 系膜钳(长头)
Mesentery forceps
Φ5X330



101.078C 系膜钳(加长中空)
Mesentery forceps
Φ5X330



101.092 穿线钳
Threading forceps
Φ10 X330



101.100 双开小抓钳
Small grasping forceps
Φ5X330



101.101 粗齿固定抓钳
Fixed grasping forceps
Φ5X330



101.101B 粗齿固定抓钳
Fixed grasping forceps
Φ5X330



101.101C 粗齿固定抓钳
Fixed grasping forceps
Φ5X330



101.102A 胃抓钳
Stomach grasping forceps
Φ5X330



101.102B 胃抓钳
Stomach grasping forceps
Φ5X330



101.102C 胃抓钳
Stomach grasping forceps
Φ5X330




101.102E 胃抓钳
Stomach grasping forceps
Φ5X330



101.103 无损伤小抓钳
Non-destructive grasping forceps
Φ5X330




101.103B 无损伤小抓钳
Atraumatic grasping forceps
Φ5X330



101.104 八爪抓钳
Octopus grasping forceps
Φ5X330



101.105 拉钩钳
hook forceps
Φ5X330




101.107 成切钳
Punch forceps
Φ5X330



101.108 无创抓钳
Atraumatic forceps
Φ5X330


其它抓钳 Other forceps (Flexible)




一次性手术器械 Disposable surgical instrument




穿刺器、转换器 Trocar, Reducer



哈巴狗抓钳 Artery/Vein Clip, Clip Applier



KE-102 微型腹腔镜手术器械 Miniature instruments for Laparoscopic surgery



KE-103 小切口手术器械 Instruments for small Incision operation

推荐标准套

产品编号	产品名称	规格	数量
CODE	NAME	SPECIFICATION	QUANTITY
103.001	镊类夹钳	Clip applicator	1
103.002	取石钳	Biopsy forceps	1
103.003	直形剪刀	Scissors(straight)	1
103.004	弯形剪刀	Scissors(curved)	1
103.005	钩形剪刀	Hook scissors	1
103.006	直形分离钳	Dissecting forceps(straight)	1
103.007	弯形分离钳	Dissecting forceps(curved)	1
103.008	直形分离钳	Dissecting forceps(90°)	1
103.009	重型抓钳	Forceps	1
103.011	L形拉钩	Hook(L-shaped)	2
103.012	S形拉钩	Hook(S-shaped)	2
103.013	活动拉钩	Movable retractor	1
103.014	冲洗管	Suction & Irrigation tube	1
KE-001B	纤维光导	Fiber optic light cable	1
103.015	电钩	Hook-electrode	1
103.017	电棒	Ballpoint-electrode	1
103.020A	推结器	Knit-pusher	1
103.021	持针钳	Needle-holder forceps	1

KE-104 胸腔镜手术器械 Instruments for thoracoscopic surgery

推荐标准套

产品编号	产品名称	规格	数量
CODE	NAME	SPECIFICATION	QUANTITY
101.001	镊类钳	Clip applicator (single)	1
101.075A	肺叶钳	Lung lobe forceps	1
104.003	宽齿固定抓钳	Fastness forceps(thick tooth)	1
104.004	宽嘴无齿抓钳	Broadness mouth forceps	1
101.067A	五叶扇形抓钳	Pantlelike fan-shaped forceps	1
104.006	无齿抓钳	Altraumatic forceps	1
104.007	长柄抓钳	Grasping forceps (long-jaw)	1
101.005	取石钳	Biopsy forceps	1
104.008	鸭嘴大抓钳	Duck-jaw forceps	1
201.021	穿刺针	Bile duct needle	1
102.020B	推结器	Knit Pusher	1
101.043A	双刃抓钳	Claw forceps	1
101.043B	双刃抓钳	Claw forceps	1
101.043D	双刃抓钳	Scissors(double-action straight)	1
101.017	弯剪刀	Scissors(double-action curved)	1
101.019	弯剪刀	Scissors(double-action curved)	1
101.022	钩形剪刀	Hook scissors	1
101.012	直形分离钳	Dissecting forceps(straight)	1
101.013	弯形分离钳	Dissecting forceps(curved)	1
101.011B	直形分离钳	Dissecting forceps(90°)	1
101.027	电钩	Hook-electrode	1
101.026	电棒	Ballpoint-electrode	1
101.028	电刀	Knife-electrode	1
101.030	电铲	Spatula-electrode	1
101.031	电针	Needle-electrode	1
104.024	穿鞘器	Trocar	1
104.025	穿鞘器	Trocar	1
104.026	穿鞘器	Trocar	1
104.027	穿鞘器	Trocar	1
104.031A	胸腔吸引管(弯头)	Irrigation tube	1
104.031B	胸腔吸引管(直头)	Irrigation tube	1
104.032	淋巴抓钳	Lymph forceps	1
101.057	抓钳	Forceps	1

KE-106 腹腔镜甲状腺手术器械 Thyroid operating instruments

推荐标准套

产品编号	产品名称	规格	数量
CODE	NAME	SPECIFICATION	QUANTITY
101.019	弯剪刀	Scissors(double-action curved)	1
101.013	弯形分离钳	Dissecting forceps(curved)	1
105.003	专用穿鞘器	Trocar	1
106.004	专用穿鞘器	Trocar	1
101.071A	转换器	Reducer	1
101.026	自动四位持针钳	Needle-holder forceps(restoration)	1
106.007	剥离器	Stripper	1
106.008	注水器	Injection needle	1
105.009	分离器	Thyroid dissecting	1
106.010	可弯剥离器	Stripper(curved)	1
106.011	专用拉钩(左、右)	Hook	2
101.068	可拆冲洗管	suction & irrigation tube	1

KE-107 腹腔镜胆道取石手术器械 Instruments of Biliary tract fetch

推荐标准套


产品编号	产品名称	规格	数量
CODE	NAME	SPECIFICATION	QUANTITY
107.001	胰管食管钳	Vein forceps	1
107.002	血管阻断钳	Vein interdiction forceps	1
107.003	胆道取石钳	Biliary tract fetch forceps	1
107.004	软管穿刺器	Trocar(hose)	1

KE-101DK单孔(经脐部)腹腔镜手术器械 Single-hole (After umbilical)


标准配套

101.094DK	穿刺通道器	φ25x40mm	1套
101.092DKa	穿刺软管	φ6x70mm	3支
101.092DKb	穿刺针	φ5.5x90mm	2支
101.093DK	简易穿刺器	φ10.5x90mm	1套
101.019DK	双动弯剪	φ5x410mm	1把
101.013DK	分离钳	φ5x410mm	2把
101.080DK	塑料钛夹钳	φ5x410mm	1把
101.023DK	持针钳	φ5x410mm	1把
101.036DK	鸭嘴抓钳	φ5x410mm	1把
101.035DK	胆囊抓钳	φ5x410mm	1把
101.068DK	吸引器	φ5x410mm	1支
101.027DK	电钩	φ5x410mm	1支
101.028DK	电棒	φ5x410mm	1支
101.095DK	多用通道器	φ5x410mm	1套
	密封帽		3只



KE-101DK单孔 (经济部) 腹腔镜手术器械 Single-hole (After umbilical)



KE-201A电动子宫切除器 Electric uterus-cutter



KE-201C 妇科腹腔镜手术器械 (选购件) Instrument of Gynaecology Laparoscope




KE-201C 妇科腹腔镜手术器械 (选购件)




KE-201C 妇科腹腔镜手术器械 (选购件) Instrument of Gynaecology Laparoscope



KE-201B 子宫切除器械 (阴式) Uterus-cutter(process vagina)




KE-202 宫腔镜 Hysteroscope




KE-203 宫腔镜(检查镜)Cavitas inspectoscope



KE-204 宫腔电切镜 Hystero Resectoscope



KE-301 尿道膀胱镜 Cysto-Urethroscope



KE-304 膀胱碎石镜 Lithotriposcope

304.002A 大力碎石钳 Strong-force lithoclast

304.001

304.002B 直咬式碎石钳 Strong-force lithoclast

304.004

304.005

301.014

301.015

推荐标准套			规格
产品编号	产品名称	规格	
CODE	NAME	SPECIFICATION	
304.001	内窥镜70°	Endoscope 70°	Φ4
304.002	大力碎石钳	Lithoclast	
301.014	进水接头	Intake connector	
301.015	出水接头	Outtake connector	
301.006	鞘套	Sheath	Φ19.5X230
304.004	冲吸器	Irrigating connection	
304.005	冲洗接头	Outtake connector	
KE-001	导光束	Fiber optic light cable	Φ4X1800

KE-305 尿道切开镜 Urethrotome Telescope

305.003

305.002

305.001

305.004

301.015

301.014

推荐标准套			规格
产品编号	产品名称	规格	
CODE	NAME	SPECIFICATION	
305.001	内窥镜12°	Endoscope 12°	Φ4x302
305.002	切开镜	Working handle	Passive
305.003	镜鞘	Sheath	
305.004	附加器	Adapter	
305.005	柳叶刀	Lance blade	
305.006	半圆刀	Hall-moon blade	
305.007	钩状刀	Sickle-shaped blade	
301.014	进水接头	Intake connector	
301.015	出水接头	Outtake connector	
KE-001A	导光束	Fiber optic light cable	Φ4x1800

KE-501 鼻窦镜 Sinoscope

501.028 镰状刀 sickle-shaped knife

501.029 镰状刀 sickle-shaped knife

501.019 鼻窦刮匙钳 Antrum curette

501.018 鼻窦刮匙钳(有孔) Antrum curette

501.027 剥离子 Detacher

501.024 直吸引管 Suction catheter, straight

501.025 直吸引管 Suction catheter, straight

501.027 弯头吸引管 Suction tube, curved

501.026 弯头吸引管 Suction tube, curved

501.031 穿刺针 Trocar

501.036 持镜架 Endoscope Handle

501.001 0° 内窥镜 Sinoscope

501.001 30° 内窥镜 Sinoscope

501.001 70° 内窥镜 Sinoscope

501.012 鼻腔组织钳(带吸引) Cutting forceps

501.004 光学钳 Optics forceps

501.016 鼻窦组织钳70° Cutting forceps

501.017 鼻窦组织钳70° Cutting forceps

501.015 鼻窦组织钳100° Cutting forceps

501.015 鼻窦组织钳 Optics forceps

501.005 直剪 Nasal scissors straight

501.006 左弯剪 Nasal scissors left

501.007 右弯剪 Nasal scissors right

501.011 鼻腔组织钳30° Polypus forceps

501.012 鼻腔组织钳45° Polypus forceps

501.013 鼻腔组织钳90° Polypus forceps

501.008 鼻咬骨钳(前咬) Cutting forceps stryken

501.009 鼻咬骨钳(后咬) Cutting forceps retrograde

501.009 鼻咬骨钳(后咬) Cutting forceps retrograde

501.027 上颌窦咬骨钳 Maxillary sinus rongeur

KE-503 鼻窦电手术刀 Nose electric-drill

501.002 手机 Debrider unit

503.004 平切刀 Filler blade

503.005 锯齿刀 Footed blade, thick

503.006 中齿刀 Footed blade, shapp


503.007 打磨刀 Burr blade

503.008 内切弯刀 Inner cutting curved blade


501.003 脚踏开关 Full-functional foot switch

501.001 控制器 Control unit

KE-502 支撑喉镜 Hold Laryngoscope




KE-506 食道镜 Esophagoscope instruments



KE-505 耳科手术器械包 Otoloscope Instruments




KE-504 小儿支气管镜 Pediatric Bronchoscopy sets



KE-401 关节镜 Arthroscope



KE-402 椎间盘镜手术器械 Endoscopic Discectomy system



KE-601 数字内镜图像显示系统 Endoscope display systems

数字内镜图像显示系统	
产品编号	产品名称
KE-601A	内镜摄像系统(基本型)
601.101	CCD彩色摄像机(单晶片)
601.102	彩色监视器松下
601.103	光学接口
601.104	稳压电源
KE-605	冷光源
KE-001A	导光束
601.105	专用仪器台车

产品规格(mm)	
单晶片	松下14"
型号可选	12V/1A
250W/氙灯	Φ4X1800
多功能架	

KE-601B 内镜摄像系统(标准型)		
产品编号	产品名称	产品规格(mm)
601.201	3CCD彩色摄像机	三晶片
601.202	彩色监视器	索尼松下14"
601.203	光学接口	型号可选
601.204	三晶片专用电源	12V/3A
KE-605	冷光源	250W/氙灯
KE-001A	导光束	Φ4X1800
601.105	专用仪器台车	多功能架

KE-601C 内镜摄像系统(高级型)		
产品编号	产品名称	产品规格(mm)
601.301	3CCD彩色摄像机	进口3CCD
601.302	彩色监视器	松下17"
601.303	专用监视仪	
601.304	三晶片光学接口	
601.204	三晶片专用电源	250W/氙灯
KE-605	冷光源	Φ4X1800
KE-001A	A导光束	Φ4X1800
601.105	专用仪器台车	多功能架

高清晰内镜图像显示器 高档3CCD内镜图像显示器

其它附属设备 Other equipment

KE-610 耳鼻喉科诊疗台系列 ENT unit

KE-611 模拟操作箱 Simulated training box

KE-603 医用加压器 (冲吸泵) Suction & Irrigation apparatus

整机示意图

■ 产品介绍 Products introduction:

- 用于腔镜手术的清洗、吸出废液及组织碎末
Washing and suction the effluent or fragment in the enterocoele surgery

■ 功能特点 Function Characteristic

- 冲洗和吸引交替循环工作
Irrigate and suction alternately
- 大流量 Large flow
- 根据需要调节冲洗和吸引压力
Adjusted irrigate and suction pressure as the need
- 防回流 Reflux prevention to prevent the effluent into the pump.

■ 技术参数 Technical Parameters

- 额定功率 Rated pwr.:50VA
- 正压力设定范围 Range of pos-pressure:10-50kPa
- 负压力设定范围 range of neg-pressure: 50--10kPa
- 噪声 Noise : ≤50dB

■ 整机示意图解

① 电源开关	⑨ 吸压力显示 (MAX)
② 冲泵出气口	⑩ 吸压力设置 (MAX)
③ 冲压力显示 (MAX)	⑪ 吸引管
④ 冲压力设置 (MAX)	⑫ 溢流阀 (内置)
⑤ 冲气管	⑬ 负压瓶
⑥ 冲液管 (可高温消毒) 2米	⑭ 瓶架
⑦ 冲气液转换针	⑮ 冲吸器 (选配)
⑧ 吸泵进气口	⑯ 冲洗液瓶 (自备)

KE-604 医用加压器 (灌注泵) Pump

整机示意图

■ 产品介绍 Products introduction:

- 用于宫腔镜的膨宫加压和泌尿外科手术的液体灌注,也可用于关节镜、腹腔镜等手术的液体灌注。
Used for Hysteroscope's swell and force,Urology surgery's liquid irrigation,also be able to use in Arthroscopy,Laparoscope surgery.

■ 功能特点 Function Characteristic

- 采用滚动挤压泵产生大流量水流,冲洗液管可以高温高压消毒,不会产生二次污染。压力可以根据需要调节设定,也可以进行实时监控。
Use roller pump to get large flow,irrigation tube can be autoclavable sterilization to avoid secondary pollution.Pressure is adjustable on request and real time monitored.

■ 技术参数 Technical Parameters

- 额定功率 Rated pwr.:100VA
- 压力设定范围 Range of press: 10-400mmHg (2-53.3kPa)
- 流量设定范围 Range of flw.:0.1-1.0L/min
- 噪声 Noise : ≤50dB

■ 适用范围
宫腔手术、尿道膀胱手术、关节手术、腹腔镜手术冲洗等

■ 整机示意图解

① 电源开关	⑨ 压力设定上调键
② 腔内实际压力显示屏	⑩ 运行键
③ 设定流量下调键	⑪ Kp, mmHg转换键
④ 设定流量显示屏	⑫ 指示灯
⑤ 设定流量上调键	⑬ 泵轮
⑥ 压力设定下调键	⑭ 出水接头
⑦ 停止键	⑮ 传感输液器
⑧ 设定压力显示屏	⑯ 容器 (用户自备)

KE-605冷光源 Medical cold light source



整机示意图

■ 产品介绍 Products introduction:

- 用于腹腔镜手术的腔内照明
Used for intracavity lighting during the enterocoelia surgery

■ 功能特点 Function Characteristic

- 亮度高,性能稳定
High luminance, steady function

■ 技术参数 Technical Parameters

- 氙灯功率 Lamp power: 350W
- 照度 Illumination: $>2000000\text{LX}$
- 色温 Color temp: $5900\text{K} \pm 10\%$
- 灯泡寿命 Lamp life: $\geq 500\text{h}$

KE-606 器械清洗成套设备 Cleaning device



整机示意图

■ 整机示意图解

- ① 电源开关
- ② 工作指示灯
- ③ 待命指示灯
- ④ 水泵按钮开关
- ⑤ 气泵按钮开关
- ⑥ 出水接口
- ⑦ 出气接口
- ⑧ 水枪枪头托盘
- ⑨ 水枪
- ⑩ 气枪
- ⑪ 带快速接头水枪软管
- ⑫ 带快速接头气枪软管

■ 产品介绍 Products introduction:

- 用于对各种手术器械、器皿以及软镜及内管道的外表面或内管道的清洗、吹干。防止手术器械二次污染。
Wash and dry various kinds of instruments, vessels' tube and surface. Avoid instruments' secondary pollution.

■ 功能特点 Function Characteristic

- 清洗压力高、震动小、噪音低、安全性能高,快速高效率。吹气气压高,储气容量大。
High pressure, little shake, low noise, safe high efficiency. High gas pressure, high storage.

■ 技术参数 Technical Parameters

- 输出功率 Output pwr.: 800VA
- 输出气压 Output gas press.: 0.4Mpa
- 输出水压 Output liquid press.: 0.5Mpa
- 噪声 Noise: $\leq 50\text{dB}$

KE-607 气腹机 Insufflator



整机示意图

■ 产品介绍 Products introduction:

- 用于腹腔镜手术中建立二氧化碳气腹,为腹腔镜手术或检查提供广阔的空间和视野。
Build CO2 space in the laparoscope surgery, for offering wide view in the surgery.

■ 功能特点 Function Characteristic

- 采用MOTORALA压力传感器,三路压力测控,大流量供气;寿命长,安全可靠。
Use MOTORALA pressure sensor, three-way pressure observe and control, large flow air supply. Long life, safe and reliable.

■ 技术参数 Technical Parameters

- 额定功率 Rated pwr.: 40VA
- 流量设定范围 Range of flw.: 10-30L
- 压力设定范围 Range of press.: 5-30mmHg
- 噪声 Noise: $\leq 50\text{dB}$

KE-608 高频电刀 High frequency electric-cutter



■ 产品介绍 Products introduction:

- POWER-420X型高频手术器是一种多功能高频电刀,广泛适用于普外、心胸、泌尿、骨科、整形外科以及配合内镜、腹腔镜手术、为闭环调幅调功型电刀。


■ 技术参数 Technical Parameters

- 仪器类型: I类CF型防除颤普通设备
- 工作频率: 512KHz
- 高频漏电流 $\leq 150\text{mA}$
- 输出功率:
纯切: 300W (500 Ω)
混切1: 250W (500 Ω)
混切2: 200W (500 Ω)
混切3: 120W (500 Ω)
电凝: 120W (500 Ω)
双极电凝: 50W (100 Ω)
- 运行凡是: 间隙加载连续运行10S/30S
- 电源: $220\text{V} \pm 22\text{V}, 50\text{Hz} \pm 1\text{Hz}$, 最大电流 $\leq 3.5\text{A}$
- 外形尺寸: $450 \times 355 \times 175$
- 净重: 11Kg, 毛重: 19kg.
- 使用条件: 环境温度 $5^\circ\text{C} - 40^\circ\text{C}$, 湿度 $\leq 80\% \text{RH}$


禁忌症:

使用心脏起搏器或有源植入物的患者; 手术过程中高频电流可能流过横截面较小的部位, 为避免不必要的凝结, 最好使用双极技术。

其它附属设备 Other equipment (Sterilising trays)



规格: 43X8.5X4.2cm



规格: 27X8.5X4.2cm



规格: 54.5X22.5X11cm

